

## ADJUSTABLE FLANGE

### NASOPHARYNGEAL AIRWAY (Latex Free)

The nasopharyngeal airway is manufactured with high quality hand-picked raw material and with smooth surface to facilitate easy insertion and removal while providing an alternative airway.

#### Features:

- Soft
- Sterile
- Disposable
- Economical
- Adjustable Flange
- Fast application

Part No.	ID (mm)	OD (mm)	Length (mm)
NS-12FR-GALFHC-DT	1.5	4.0	109
NS-14FR-GALFHC-DT	2.0	4.7	109
NS-16FR-GALFHC-DT	2.5	5.3	109
NS-18FR-GALFHC-DT	3.2	6.0	109
NS-20FR-GALFHC-DT	3.9	6.7	124
NS-22FR-GALFHC-DT	4.8	7.3	136
NS-24FR-GALFHC-DT	5.5	8.0	169
NS-26FR-GALFHC-DT	6.3	8.7	169
NS-28FR-GALFHC-DT	6.8	9.3	169
NS-30FR-GALFHC-DT	7.0	10.0	169
NS-32FR-GALFHC-DT	7.6	10.7	169
NS-34FR-GALFHC-DT	8.6	11.3	169
NS-36FR-GALFHC-DT	9.0	12.0	169

## ROBERTAZZI DESIGN

### NASOPHARYNGEAL AIRWAY (Latex Free)

The Robertazzi Design Nasopharyngeal Airway is designed and manufactured to be soft and flexible to minimize discomfort. Only quality neoprene material is used. Comes with smoothly rounded edges to provide more comfort to patients.

#### Features:

- Soft
- Sterile
- Disposable
- Economical
- Trumpet Flange
- Fast Application

Part No.	ID (mm)	OD (mm)	Length (mm)
NS-20FR-GRLFHC-DT	3.9	6.7	114
NS-22FR-GRLFHC-DT	4.8	7.3	114
NS-24FR-GRLFHC-DT	5.5	8.0	121
NS-26FR-GRLFHC-DT	6.3	8.7	127
NS-28FR-GRLFHC-DT	6.8	9.3	133
NS-30FR-GRLFHC-DT	7.0	10.0	140
NS-32FR-GRLFHC-DT	7.6	10.7	146
NS-34FR-GRLFHC-DT	8.6	11.3	159
NS-36FR-GRLFHC-DT	9.0	12.0	171

