



## MedEvac4 Backpack



**CE** This product conforms to the Medical Device Regulations 2017/745

The MedEvac4 Backpack is designed to carry our MedEvac4 Tactical Stretcher as well as medical supplies and first aid equipment used by the combat medic. It is ergonomically designed for fast access with slots for labelling. The back of the stretcher is outfitted with 2 metal strips to prevent the stretcher from knocking against the spine of the combat medic when running.

### Key features at glance:

- Carries 1 MedEvac4 Tactical Stretcher
- Lightweight
- Transparent pockets for labelling of supplies
- Rugged and Durable
- Several pockets for storing medical supplies
- Metallic frame backing protects the medic's back

Specifications:	
Dimension:	62cm x 42cm x 25cm
Net weight:	1.2kg
Accessories:	MedEva4 Tactical Stretcher Tactical IV Pole Medical Supplies

\*Specifications are subject to change without prior notice.



[www.MedAid.ie](http://www.MedAid.ie)

### Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478  
Email: [info@goldenseason.com.sg](mailto:info@goldenseason.com.sg) Website: [www.goldenseason.com.sg](http://www.goldenseason.com.sg) | [www.gseason.com](http://www.gseason.com) Co. Regn No: 03734/1984E  
GST Registration No. M2-0067495-2