



## Stretcher Wheel Cart/Wheeled Litter Carrier

(PATENTED)

**CE** This product conforms to the Medical Device Regulations 2017/745



Golden Season's Stretcher Wheel Cart/Wheeled Litter Carrier is a lightweight casualty transport system that allows one person to transport a casualty strapped to a stretcher and also performs as a litter stand in the battlefield or medical stations.

The revolutionary design is foldable, and can also be totally disassembled for compact storage, easy shipping and transportation. The Stretcher Wheel Cart/Wheeled Litter Carrier is compatible with modern combat litters and is equipped with non-puncture wheels. A single user can operate the Stretcher Wheel Cart/Wheeled Litter Carrier.



[www.MedAid.ie](http://www.MedAid.ie)

### Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478  
Email: [info@goldenseason.com.sg](mailto:info@goldenseason.com.sg) Website: [www.goldenseason.com.sg](http://www.goldenseason.com.sg) | [www.gseason.com](http://www.gseason.com) Co. Regn No: 03734/1984E  
GST Registration No. M2-0067495-2



Technical Information	
Deployed dimensions:	L143 x W70 x H92 (cm)
Half Folded dimensions:	L126 x W70 x H50 (cm)
Fully Folded dimensions:	L107 x W70 x H25 (cm)
Weight:	17.5kgs excluding litter

\*Specifications are subject to change without prior notice.



## Characteristics

- ✓ Foldable and compactable
- ✓ Able to perform as a litter stand
- ✓ Lightweight – 17.5kgs only
- ✓ Pneumatic non puncture wheels
- ✓ Can be fully disassembled for easy shipping

The Stretcher Wheel Cart/Wheeled Litter Carrier is compatible with modern stretcher systems including the Medevac4, TalonII 90C, Raven 90C, T3 Assault Litters and other stretchers.

## Golden Season Pte Ltd