



Stretcher Back Rest



CE This product conforms to the Medical Device Regulations 2017/745

Golden Season's MedEvac4 Stretcher Backrest is designed for MedEvac4 Tactical Stretcher that conforms to the patient body contour to relieve discomfort and help the patient to eat while lying on the MedEvac4 Tactical Stretcher. The MedEvac4 Stretcher Backrest is also built from Altrunium.

Our MedEvac4 Stretcher Backrest is weighing approximately 1.27 kg. Easy to install to MedEvac4 Tactical Stretcher.

Specifications	
Length:	80 cm
Width:	87 cm
Height:	3 cm
Weight:	1.27 kg
Material Type:	Lightweight Aluminium
Accessories:	MedEvac4 Litter Stand, MedEvac4 Folding Stretcher

*Specifications are subject to change without prior notice.



www.MedAid.ie

Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478
Email: info@goldenseason.com.sg Website: www.goldenseason.com.sg | www.gseason.com Co. Regn No: 03734/1984E
GST Registration No. M2-0067495-2