



Field Hospital Bed Carrier



Field Hospital Bed Carrier is designed to carry the Field Hospital Bed from one place to another location. The advantages of Field hospital bed carrier are having 4 wheels on its bottom with two hand grip on both sides, this will help the medical staff for transporting the field hospital bed without lifting.

Key features at glance

- Rugged Design
- Heavy Duty Wheels
- Strong Hand Grip
- 6 straps to secured the field hospital bed.
- Inside pocket for Meal tray and Extendable IV Pole and other accessories.

Specifications	
Dimension:	103cm x 67cm x 28cm
Nett Weight:	3.11kg
Gross Weight:	27.96kg (Including Field Hospital Bed, Meal tray and IV Pole)
Hand Grip:	2 Hand Grips
Straps:	6 straps to secured the field hospital bed.

*Specifications are subject to change without prior notice.



www.MedAid.ie

Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478
Email: info@goldenseason.com.sg Website: www.goldenseason.com.sg | www.gseason.com Co. Regn No: 03734/1984E
GST Registration No. M2-0067495-2