



Tactical IV Pole



CE This product conforms to the Medical Device Regulations 2017/745

The “Extendable” IV Pole can be easily inserted into any of the 6 slots of the MedEvac4 Tactical Stretcher as well as compatible models. It serves as a stable drip stand.

Key features at glance:

- Lightweight
- Adjustable Length
- Acts as a IV Stand on a combat Stretchers
- Simple twist to lock

Specifications:	
Dimension (Retracted):	64cm x 22cm
Dimension (Extended):	110cm x 22cm
Pole diameter:	1.2cm
Weight:	200g
Accessories:	MedEvac4 Tactical Stretcher MedEvac4 Litter Stand MedEvac4 Carrier / Backpack

*Specifications are subject to change without prior notice.



www.MedAid.ie

Golden Season Pte Ltd