



MedEvac4 STANAG 2040 NATO Stretcher



The MedEvac4 Fold STANAG 2040 Stretcher is the NATO version of the reliable and proven MedEvac4 Stretcher Platform and is manufactured in accordance with NATO Standards (STANAG 2040). It folds into 4 for compact transport / storage, and yet maintains the extremely high load bearing capability that our stretchers are renowned for.

CE This product conforms to the Medical Device Regulations 2017/745

Specifications:	
Dimensions Extended	2290mm X 584mm X 171mm
Dimensions Folded	630mm X 215mm X 220mm
Canvas Specification	Warp: Minimum 2400N / 50mm Weft: Minimum 1400N / 50mm Compliant with MIL-L-16446A Standards
Frame	Aerospace Grade Aluminum
Handles	Retractable, Automatic locking, Ergonomically Designed
Weight	7kg
CE Mark	Yes

*Specifications are subject to change without prior notice.



www.MedAid.ie

Golden Season Pte Ltd