



Side Rails for Mark II ICU Bed



UPSB Side Rail is designed to fit in for Mark II Field Hospital Bed. It is designed to secure the patient from falling down on the field hospital bed and serve as arm rest for the patient to prevent discomfort while lying on the field hospital bed.

Key features at glance

- Manufactured using Altrunium®
- Heavy duty side rail
- Design to fit in for Mark II Field Hospital Bed
- Automatic locking mechanism
- Easy to use.

Specifications	
Length:	74.5cm
Width:	3cm
Height:	30cm
Material:	Built from Altrunium®
Accessories:	Mark II Field Hospital Bed (Sold separately)

*Specifications are subject to change without prior notice.



www.MedAid.ie

Golden Season Pte Ltd