



## Single Fold STANAG 2040 NATO Stretcher



Golden Season's Single Fold STANAG 2040 Stretcher serves as a reliable casualty evacuation platform and is manufactured in accordance with NATO Standards (STANAG 2040).

**CE** This product conforms to the Medical Device Regulations 2017/745

Specifications:	
Dimensions Extended	2290mm X 570mm X 153mm
Dimensions Folded (Handles Retracted)	1980mm X 85mm X 153mm
Canvas Specification	Warp: Minimum 2400N / 50mm Weft: Minimum 1400N / 50mm Compliant with MIL-L-16446A Standards
Frame	Aerospace Grade Aluminum
Handles	Retractable, Automatic locking, Ergonomically Designed
Weight	7.5kg
CE Mark	Yes

\*Specifications are subject to change without prior notice.



[www.MedAid.ie](http://www.MedAid.ie)

### Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478  
Email: [info@goldenseason.com.sg](mailto:info@goldenseason.com.sg) Website: [www.goldenseason.com.sg](http://www.goldenseason.com.sg) | [www.gseason.com](http://www.gseason.com) Co. Regn No: 03734/1984E  
GST Registration No. M2-0067495-2